**GB STÉRIMAR™**  
**HYPERTONIC SEA WATER SOLUTION MICRODIFFUSION**

**The nose is the key to a good breathing.**

The nose and the nasal mucus play an important role. They humidify and warm respired air acting as an air conditioner. When it is clogged, the nose is no longer able to properly carry out its regulatory functions. In addition, the mucus can often be a source of contamination, which can result in ENT disorders because it is a carrier of pathogenic germs.

**Since 1975, STÉRIMAR™, pioneer of sea water nasal spray, has put its expertise of the sea and its natural resources at the service of nasal health.**

STÉRIMAR™ is medical device used and suitable for the whole family from babies to adults. Filtered, 100% natural and preservative-free, STÉRIMAR™ sea water benefits from the recognised properties of mineral salts and marine trace elements (zinc, copper, manganese, iron, magnesium, selenium...)

STÉRIMAR™ thanks to the fineness of its microdiffusion, allows you to gently wash the nasal cavities and facilitates the absorption of trace elements whilst avoiding the spread of mucus.

STÉRIMAR™ respects the environment: this aluminum container can be recycled. Its packaging is made of recycled cardboard and its package insert is printed on environmentally friendly paper without chlorine. The sea water is drawn from the Atlantic coast and after checking the physiological purity, it is made for use.

**DIRECTIONS:**

STÉRIMAR™ can be used in ALL POSITIONS and with maximum security. The anatomical applicator is suitable for small infant nostrils as well as the nostrils of an adult.

**Careful.** Pressurised container: May burst if heated. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not pierce or burn, even after use. Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F. Keep out of reach of children.

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- Decongests
- Helps to fight colds

### INDICATION OF USE:

- **Quickly decongest** the nose,
- **Eliminate mucus**, even thick,
- **Help to fight colds and limit the risks of secondary infection** by washing away nasal cavities. Nasal cleansing (if blocked nose) improves the effectiveness of local medicated treatments. **If the nasal mucosa is irritated, feelings of tingling may occur when using.**

### RECOMMENDED USE:

For babies from 3 months to 3 years. *On first use, spray once to start the system.*

Gently insert the nozzle into the nostril. Spray briefly 1 to 2 times in each nostril. Let any excess solution run out before wiping the nose.

**For a deeper cleansing of the nasal cavities:** Place the baby from 0-2 years old lying on his/her back or sitting up from the age of 2. Baby's head tilted to the side (never backwards), gently insert the nozzle into the uppermost nostril. Press the nozzle for approximately 2 seconds. Let any excess solution run out of the other nostril then wipe the nose. Repeat the operation in the other nostril if necessary.

After each use clean the nozzle with soapy water and dry.

**Use 2 to 6 times a day or more until decongestion.** Intended for temporary use if blocked nose.

### COMPOSITION:

Sea water: 60 ml (salt: approx. 19g/l)

Purified water qsp 100ml

Pentahydrate copper salt

Monohydrate manganese salt

Net volume: **50ml** e  **89**

Approximately 150 sprays

Aseptically packaged.

Pressurised container with an inert gas (nitrogen).

The solution is held within an inner container and does not come into contact with the nitrogen.

\*Solution with a higher salt concentration than the human body's cells.

Sea water HYPERTONIC\*

Salts between 17.5 and 20.5 g/l

Enriched with COPPER and MANGANESE.

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